DE B. Braun Avitum AG
Schwarzenberger Weg 73-79
34212 Melsungen - Germany

BE B. Braun Medical N.V./S.A.
1831 Diegem - Belgium

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CS B. Braun Medical, s.r.o.
148 00 Praha 4 - Czech Republic

SK B. Braun Medical s.r.o.
831 03 Bratislava - Slovakia

GB B. Braun Medical Ltd.
Sheffield, S35 2PW
United Kingdom

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


* For each country a specific letter is added at the end of the reference

**Available only in Spain

EN	Do not reuse
DE	Nicht wiederverwenden
ES	No reutilizar
FR	Ne pas réutiliser
IT	Non riutilizzare
NL	Niet hergebruiken
CS	Nepoužívejte opakovaně
DA	Må ikke genanvendes
MY	Jangan gunakan semula
PT	Não reutilizar
RO	A nu se refolosi
RU	Повторное использование запрещено
SK	Nepoužívajte opakovane


EN	Consult instructions for use
DE	Gebrauchsanweisung beachten
ES	Consúltese las instrucciones de uso
FR	Se reporter à la notice d'utilisation
IT	Consultare le istruzioni per l'uso
NL	Raadpleeg de gebruiksaanwijzing
CS	Přečtěte si návod k použití
DA	Se brugsanvisningen
MY	Rujuk arahan penggunaan
PT	Consulte as instruções de utilização
RO	Consultați instrucțiunile de utilizare
RU	См. инструкцию по применению
SK	Pozrite návod na použitie

EN	Manufacturer
DE	Hersteller
ES	Fabricante
FR	Fabricant
IT	Fabbricante
NL	Fabrikant
CS	Výrobce
DA	Fabrikant
MY	Pengilang
PT	Fabricante
RO	Producător
RU	Изготовитель
SK	Výrobca

EN	Date of manufacture
DE	Herstellungsdatum
ES	Fecha de fabricación
FR	Date de fabrication
IT	Data di fabbricazione
NL	Productiedatum
CS	Datum výroby
DA	Fremstillingsdato
MY	Tarikh pengilangan
PT	Data de fabrico
RO	Data fabricației
RU	Дата изготовления
SK	Dátum výroby

EN Instructions for use

Carefully read the Instructions for use before usage

Device description

Actreen® Lite Cath is a pre-lubricated, sterile and ready to use catheter.

Composition

Actreen® Lite Cath is made of PVC (catheter), Ethylene Vinyl Acetate (EVA) (prehension part), UV acrylate or cyanoacrylate (catheter and the prehension part), the lubricant is water and glycerin based and the primary packaging is made of Polyethylene (PE).

Not made with DEHP and Latex.

Sterility

Sterilized using Ethylene Oxide.

Intended purpose

The intended use (purpose) of Actreen® Lite Cath is to provide intermittent urinary catheterisation of the bladder, by introducing the catheter through the urethra to drain the urine from the bladder.

Indication

Actreen® Lite Cath is indicated for patients who need intermittent urinary catheterisation due to temporary or chronic urine retention or voiding dysfunction.

Patient population

Actreen® Lite Cath can be used in male and female adult or paediatric patients with a body weight superior or equal to 3,5 kg.

The responsible healthcare professional shall prescribe the appropriate size, tip and length of the catheter as well as the frequency of intermittent catheterisation based on patient's age, gender, characteristics, and conditions (body profile in general, pregnancy, paediatric population) according to state-of-the-art medical knowledge.

For pregnant women, the responsible healthcare professional shall assess the appropriate length and diameter of the catheters, since in pregnancy the length of the urethra is altered.

For paediatric patients, the capability to perform intermittent self-catheterisation should be determined by the responsible healthcare professional according to the state-of-the-art medical knowledge.

Intended User

The intended users are healthcare professionals, caregivers, and patients.

Caution: The responsible healthcare professional must ensure that the patients and caregivers are duly trained and instructed in intermittent (self) catheterisation with the catheter, as use of the device by an untrained user may cause trauma or compromise the sterility of the product, potentially leading to infection.

Intended Use Environment

Home, daily life activities environments, hospitals and healthcare facilities.

Contraindications

Contraindication related to the product:

Known hypersensitivity reaction to any material used (see composition).

Absolute contraindications to intermittent catheterisation:

High intravesical pressure that would require continuous free drainage to avoid renal damage. Urinary fistula. Mitrofanoff procedure. Priapism. After urethral surgery. Suspected urethral tumor. Pelvic / urethral injury.

Relative contraindications to intermittent catheterisation with special cautions:

Please contact your healthcare professional before using this product if any of the following conditions are present: Abnormal urethral anatomy (e.g., stricture), false passage, bladder neck obstruction, diseases of penis (e.g., injury, tumor, or infections), after gender reassignment surgery, tendency for autonomic dysreflexia with bladder filling/urethral instrumentation (despite treatment), urethral bleeding, prostatic or bladder neck surgery, female genital mutilation, stent or artificial prosthesis in the urogenital tract, bladder capacity < 200 ml, or high fluid intake regimen.

Contraindication for intermittent self-catheterisation:

Inability to self-catheterize in the absence of an appropriately trained caregiver or attendant for patient with poor manual dexterity and cognitive impairment.

Warning: If you have any of the above-named contraindications or are in doubt, please contact your responsible physician immediately before the next catheterisation on how to proceed.

The responsible physician must define contraindications for an intermittent catheterisation based on patient characteristics and conditions according to state-of-the-art medical knowledge.

Residual Risks

Residual risks are included in other chapters (e.g., warnings etc.).

Side effects

Side effects related to the product

Unknown hypersensitivity reaction to any material used leading to anaphylactic or anaphylactoid reactions.

Side effects related to intermittent catheterisation

Catheter-associated Urinary Tract Infection (CAUTI), in some instance associated with pyelonephritis.

Urethritis / epididymo-orchitis/ prostatitis.

Trauma of the urinary tract including urethral trauma / urethral bleeding / false passage / urethral strictures (predominantly in men) / meatal erosion and stenosis / bladder perforation.

Discomfort, soreness, pain.

Bladder calculus, stone formation.

Warnings: If severe hypersensitivity reactions occur, remove the product immediately, discontinue its use and contact a physician.

If you experience reduced or interrupted urine flow or that the catheter does not drain the bladder properly, contact a physician as it may indicate that you should use another size or type of catheter.

If you have symptoms of CAUTI or infection of urogenital tract: fever, chills, pyuria, dysuria, urgency, cloudy urine, suprapubic pain, discomfort in emptying the bladder, blood in the urine, increased bladder spasms, autonomic dysreflexia, mucus, contact a healthcare professional.

If you have any signs of trauma, discomfort when emptying the bladder, bleeding, a frequent urge to urinate or blood in the urine, discontinue use immediately and contact a healthcare professional.

Ungentle insertion and removal can cause trauma and aggravate discomfort and pain. If you have doubts about the catheter handling, contact a healthcare professional.

Avoid introduction of pubic hair. Bladder calculus or stone formation are related to the introduction of pubic hair and presence of mucus.

Cautions

Fluid intake should be adequate unless contraindicated. Patient should avoid excessive fluid intake that may cause bladder overdistension.

Constipation may cause pressure preventing drainage from catheter: Maintain regular bowel function.

For all intended users, a trained period is necessary to be autonomous on catheterisation.

The healthcare professional must ensure that the patients and caregivers are duly trained and instructed in intermittent (self) catheterisation with the catheter, as use of the device by an untrained user may cause trauma or compromise the sterility of the product, potentially leading to infection.

Caution specific for healthcare professionals: Apply with caution if the patient produces urine with many particles clearly distinguishable by the naked eye, as it may lead to transient urine retention.

Catheterisation should take place as quickly as possible after opening the package as it may promote contamination or compromise product characteristics (e.g., drying lubricant).

Do not fold the catheter as it may potentially promote trauma to the urinary tract.

Avoid contact between the urinary meatus and the connector or the sheath as it may potentially promote trauma to the urinary tract.

Warning

Any contamination of the catheter needs to be avoided by proper handling.

Do not use a female size catheter for an adult male patient.

Do not use a male Tieman catheter for female patient.

Cleaning, disinfection, resterilization may damage the product and creates a potential risk to the patient. It may lead to impairment of functionality or contamination of the device, potentially leading to infection.

Do not use the catheter if you observe any damage (e.g., kinked) to, blockage of, or foreign bodies in the device.

Do not use if side eyes are missing as the device may not perform as intended.

For urethral use only.

In case of suspected catheters knotting, please contact your healthcare professional immediately.

Do not use if package is damaged or after expiry date, as sterility of the product may be compromised, potentially leading to infection.

In rare cases, infections from intermittent catheterisations can lead to serious illness or even death. To reduce this risk, it is essential to follow proper handling procedures and maintain strict hygiene practices. Impairment of functionality of the device may lead to additional risks of injury, or illness of the patient.

Notice: There are additional, specific warnings in other sections of this Instruction for Use. These specific warnings must be considered as well.

Operating Instructions

1. PREPARATION:

- Carefully wash hands before and after catheterisation and clean the area surrounding the urethral opening.

If a caregiver is assisting catheterisation, local guidelines must be applied.

- Before opening, check the product is intact (no packaging damage, no visual defects, kink on tube...) Make sure the lubricant is evenly spread over the catheter. If not, spread the lubricant by touching the catheter over the packaging to allow an even distribution.

Unevenly spread lubricant may potentially promote trauma of the urinary tract.

- During opening, keep the product horizontally and remove protective film slowly to avoid lubricant loss.

2. PRODUCT OPENING:

For safety and hygienic reasons, we advise you to use the product in shortest possible time after opening.

Several techniques are available/could be applied to remove the protective film.

Caution: Whatever the technique employed pay attention not to put the catheter tube in contact with hand or objects/surfaces.

Warning: Wrong manipulation of the catheter could lead to contamination and/or impairment of functional capability of the device.

NO TOUCH TECHNIQUES

FOR THE FOLDABLE NO TOUCH TECHNIQUE

- Option 1: Grip and tear off the protective film at mark «2» to remove film from the connector (1A) and then at mark "A" to remove film from the catheter (2A).
- Option 2: Hold the catheter in one hand and with the other hand insert your thumb in the eyelet at mark «1» to remove the protective film on the connector (1B) and then hold the catheter in one hand and with the other hand insert your thumb in the eyelet at mark «A» to remove the protective film on the catheter (2B).
- For the 2 options, then fold down the sheath into accordion folds without touching the catheter (3).

Option 2 may be easier if you have reduced dexterity.

For the male length catheter, if you want to release more of the catheter open the product at mark "B" instead of mark "A".

IF YOU WANT TO RELEASE THE WHOLE CATHETER:

- After removing the film at mark «2» from the connector side, tear off the protective film at mark "C" for male length catheter and mark "B" for female length catheters, hold the connector and lift the catheter out of the packaging.

3. CATHETERISATION:

Actreen® Lite Cath has a conical staged connector connectable to the urine collection bag Urimed Bags from B. Braun.

Warning: bad connexion could lead to urine leak or disconnexion during catheterisation.

Precaution: Pay attention not to fold the pouch/bag to allow a correct flow of urine inside.

MALE (4)

- With one hand hold the penis upwards to reduce curvature of the urethra.
- With the other hand carefully introduce the catheter into the urethra until it reaches the bladder. Point your penis downwards in case of resistance. Continue inserting the catheter until urine starts to flow.
- Tiemann/Coudé catheter : In order to correctly guide the catheter, ensure that the mark on the connector is pointing in the required direction.

Warning: The curved tip of the catheter shall point upwards. Wrong orientation of the tip may lead to trauma of the urinary tract.

FEMALE (5)

- With one hand spread the outer lips of the vulva to locate the urethral opening. With the other hand carefully introduce the catheter into the urethra until it reaches the bladder and the urine drains.

MALE / FEMALE

- Manual pressure on the lower abdomen towards the end of catheterisation will ensure that the bladder has been completely emptied.

- Immediately after the end of bladder emptying, withdraw the catheter.
- To stop the catheter from dripping, turn up the connector end of the catheter when removing.

Duration of use

For single use only.

Warning: Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination of the device may lead to infection.

Catheters must be removed immediately after completion of urinary drainage.

Disposal

Discard the catheter according to the current local regulations.

Storage and handling conditions

Keep away from sunlight and heat. Keep dry. The product must be stored flat.

Warning: Do not store the product (unit) in contact with sharp materials.

Warning: The expiry date only applies when Actreen® Lite Cath is correctly stored in an undamaged packaging.

Do not use after the expiry date.

Clinical benefits & Performances

Actreen® Lite Cath allows the drainage of urine from the bladder through the urethra.

Notice to the user

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

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