Actreen® Hi-Lite Cath



EN Instructions for use
DE Gebrauchsinformation
ES Instrucciones de uso
FR Mode d'emploi
IT Istruzioni per l'uso
NL Gebruiksaanwijzing
AR تعليمات الاستخدام
CS Návod k použití

CN 使用说明书
DA Brugsanvisning
EL Οδηγίες χρήσης
FI Käyttöohjeet
LV Lietošanas norādījumi
NO Bruksanvisning
PL Instrukcja użycia
PT Instruções de utilização

RO Instrucțiuni de utilizare
RU Инструкция по применению
SK Návod na použitie
SV Bruksanvisning

Product range

			CHARRIÈRE (Ø MM)	REFERENCES
Actreen® Hi Lite Cath Nelaton 20 cm			CH06 (2,0 mm)	238306*
			CH08 (2,7 mm)	238308*
			CH10 (3,3 mm)	238310*
			CH12 (4,0 mm)	238312*
			CH14 (4,7 mm)	238314*
			CH16 (5,3 mm)	238316*
Actreen® Hi Lite Cath Nelaton 41 cm			CH08 (2,7 mm)	238208*
			CH10 (3,3 mm)	238210*
			CH12 (4,0 mm)	238212*
			CH14 (4,7 mm)	238214*
			CH16 (5,3 mm)	238216*
			CH18 (6,0 mm)	238218*
Actreen® Hi Lite Cath Tiemann 41 cm			CH08 (2,7 mm)	238108*
			CH10 (3,3 mm)	238110*
			CH12 (4,0 mm)	238112*
			CH14 (4,7 mm)	238114*
			CH16 (5,3 mm)	238116*
			CH18 (6,0 mm)	238118*

* For each country a specific letter is added at the end of the reference

EN Instructions for use

Carefully read the Instructions for use before usage

Device description

Actreen® Hi-Lite Cath is a pre-lubricated, sterile, and ready to use catheter.

Composition

Actreen® Hi-Lite Cath is made of Thermoplastic Polyolefin (TPO) (catheter), Ethylene Vinyl Acetate (EVA) (prehension part), UV acrylate or cyanoacrylate (catheter and prehension part connection), the lubricant is water and glycerin based, and the primary packaging is made of Polyethylene (PE) and Polypropylene (PP). Not made with PVC, DEHP, or Latex.

Sterility

Sterilized using Beta Irradiation.

Intended purpose

The intended use (purpose) of Actreen® Hi-Lite Cath is to provide intermittent urinary catheterisation of the bladder, by introducing the catheter through the urethra to drain the urine from the bladder.

Indication

Actreen® Hi-Lite Cath is indicated for patients who need intermittent urinary catheterisation due to temporary or chronic urine retention or voiding dysfunction.

Patient population

Actreen® Hi-Lite Cath can be used in male and female adult or paediatric patients with a body weight superior or equal to 3,5 kg.

The responsible healthcare professional shall prescribe the appropriate size, tip and length of the catheter as well as the frequency of intermittent catheterisation based on patient's age, gender, characteristics, and conditions (body profile in general, pregnancy, paediatric population) according to state-of-the-art medical knowledge. For pregnant women, the responsible healthcare professional shall assess the appropriate length and diameter of the catheters, since in pregnancy the length of the urethra is altered.

For paediatric patients, the capability to perform intermittent self-catheterisation should be determined by the responsible healthcare professional according to the state-of-the-art medical knowledge.

Intended User

The intended users are healthcare professionals, caregivers, and patients.

Caution: The responsible healthcare professional must ensure that the patients and caregivers are duly trained and instructed in intermittent (self) catheterisation with the

catheter, as use of the device by an untrained user may cause trauma or compromise the sterility of the product, potentially leading to infection.

Intended Use Environment

Home, daily life activities environments, hospitals and healthcare facilities.

Contraindications

Contraindication related to the product:

Known hypersensitivity reaction to any material used (see composition).

Absolute contraindications to intermittent catheterisation:

High intravesical pressure that would require continuous free drainage to avoid renal damage. Urinary fistula. Mitrofanoff procedure. Priapism. After urethral surgery. Suspected urethral tumor. Pelvic / urethral injury.

Relative contraindications to intermittent catheterisation with special cautions:

Please contact your healthcare professional before using this product if any of the following conditions are present: Abnormal urethral anatomy (e.g., stricture), false passage, bladder neck obstruction, diseases of penis (e.g., injury, tumor, or infections), after gender reassignment surgery, tendency for autonomic dysreflexia with bladder filling/urethral instrumentation (despite treatment), urethral bleeding, prostatic or bladder neck surgery, female genital mutilation, stent or artificial prosthesis in the urogenital tract, bladder capacity < 200 ml, or high fluid intake regimen.

Contraindication for intermittent self-catheterisation:

Inability to self-catheterize in the absence of an appropriately trained caregiver or attendant for patient with poor manual dexterity and cognitive impairment.

Warning: If you have any of the above-named contraindications or are in doubt, please contact your responsible physician immediately before the next catheterisation on how to proceed. The responsible physician must define contraindications for an intermittent catheterisation based on patient characteristics and conditions according to state-of-the-art medical knowledge.

Residual Risks

Residual risks are included in other chapters (e.g., warnings etc.).

Side effects

Side effects related to the product

Unknown hypersensitivity reaction to any material used leading to anaphylactic or anaphylactoid reactions.

Side effects related to intermittent catheterisation

Catheter-associated Urinary Tract Infection (CAUTI), in some instance associated with pyelonephritis. Urethritis/ epididymo-orchitis/ prostatitis. Trauma of the urinary tract including urethral trauma / urethral bleeding / false passage / urethral strictures

(predominantly in men) / meatal erosion and stenosis / bladder perforation. Discomfort, soreness, pain. Bladder calculus, stone formation.

Warning: If severe hypersensitivity reactions occur, remove the product immediately, discontinue its use and contact a physician.

If you experience reduced or interrupted urine flow or that the catheter does not drain the bladder properly, contact a healthcare professional as it may indicate that you should use another size or type of catheter.

If you have symptoms of CAUTI or infection of urogenital tract: fever, chills, pyuria, dysuria, urgency, cloudy urine, suprapubic pain, discomfort in emptying the bladder, blood in the urine, increased bladder spasms, autonomic dysreflexia, or mucus, contact a healthcare professional.

If you have any signs of trauma, discomfort when emptying the bladder, bleeding, a frequent urge to urinate or blood in the urine, discontinue use immediately and contact a healthcare professional.

Ungentle insertion and removal can cause trauma and aggravate discomfort and pain. If you have doubts about the catheter handling, contact a healthcare professional.

Avoid introduction of pubic hair. Bladder calculus or stone formation are related to the introduction of pubic hair and presence of mucus.

Cautions

Fluid intake should be adequate unless contraindicated. Patient should avoid excessive fluid intake that may cause bladder overdistension. Constipation may cause pressure preventing drainage from catheter: Maintain regular bowel function. For all intended users, a training period is necessary to be autonomous on catheterisation. The healthcare professional must ensure that the patients and caregivers are duly trained and instructed in intermittent (self) catheterisation with the catheter, as use of the device by an untrained user may cause trauma or compromise the sterility of the product, potentially leading to infection. Caution specific for healthcare professionals: Apply with caution if the patient produces urine with many particles clearly distinguishable by the naked eye, as it may lead to transient urine retention. Catheterisation should take place as quickly as possible after opening the package as it may promote contamination or compromise product characteristics (e.g., drying lubricant). Do not fold a non-foldable catheter (short hi-lite and CH18) as it may potentially promote trauma to the urinary tract. Avoid contact between the urinary meatus and the connector or the sheath as it may potentially promote trauma to the urinary tract.

Warning

Any contamination of the catheter needs to be avoided by proper handling. Do not use a female size catheter for an adult male patient. Do not use a male Tieman catheter for a female patient. Cleaning, disinfection, resterilization may damage the product and creates a potential risk to the patient. It may lead to impairment of functionality or

contamination of the device, potentially leading to infection. Do not use the catheter if you observe any damage (e.g., kinked) to, blockage of, or foreign bodies in the device. Do not use if side eyes are missing as the device may not perform as intended. For urethral use only. In case of suspected catheters knotting, please contact your healthcare professional immediately. Do not use if package is damaged or after expiry date, as sterility of the product may be compromised, potentially leading to infection.

In rare cases, infections from intermittent catheterisations can lead to serious illness or even death. To reduce this risk, it is essential to follow proper handling procedures and maintain strict hygiene practices.

Impairment of functionality of the device may lead to additional risks of injury, or illness of the patient.

Notice: There are additional, specific warnings in other sections of this Instruction for use. These specific warnings must be considered as well.

Operating Instructions

1. PREPARATION:

- Carefully wash hands before and after catheterisation and clean the area surrounding the urethral opening.
 - If a caregiver is assisting catheterisation, local guidelines must be applied.
 - Before opening, unfold the catheter (if applicable) and check if the product is intact (no packaging damage, no visual defects, kink on tube...) Make sure the lubricant is evenly spread over the catheter.
- If not, spread the lubricant by touching the catheter over the packaging to allow an even distribution.

Unevenly spread lubricant may potentially promote trauma of the urinary tract.

- During opening, keep the product horizontally and remove protective film slowly to avoid lubricant loss.

2. PRODUCT OPENING:

For safety and hygienic reasons, we advise you to use the product in shortest possible time after opening. Several techniques are available/could be applied to remove the protective film.

Caution: Whatever the technique employed pay attention not to put the catheter tube in contact with hand or objects/surfaces.

Warning: Wrong manipulation of the catheter could lead to contamination and/or impairment of functional capability of the device.

No touch techniques

Remove the protective film on the connector side at the mark «1».(A)

- For the foldable no touch technique:




Option 1: Grip and tear off the protective film at mark «A».(B)

Option 1

Option 2

Option 2: Hold the catheter in one hand and with the other hand insert your thumb in the eyelet at mark «A» to remove the protective film*(C)

For both options, then fold down the sheath into accordion folds without touching the catheter. Option 2 may be easier if you have reduced dexterity.(D)



• For a simplified no touch technique:

Release a movable no touch sheath by gripping and tearing off the film at marks «1» and «B». This movable part will protect the catheter when handling.(E)

If you want to release the whole catheter:

Remove the protective film on the catheter side.(F)

After removing the film at mark «1» from the connector side, hold the cone and lift the catheter out of the packaging.(G)


3. CATHETERISATION:

Actreen® Hi-Lite Cath has a conical staged connector that is connectable to the urine collection bag Urimed Bags from B. Braun.

Warning: Bad connection could lead to urine leak or disconnection during catheterisation.

Precaution: Pay attention not to fold the pouch/bag to allow a correct flow of urine inside.

MALE

- With one hand hold the penis upwards to reduce curvature of the urethra.
- With the other hand carefully introduce the catheter into the urethra until it reaches the bladder. Point your penis downwards in case of resistance. Continue inserting the catheter until urine starts to flow.
- Tiemann/Coude catheter: In order to correctly guide the catheter, ensure

that the mark on the connector is pointing in the required direction.

Warning: The curved tip of the catheter shall point upwards. Wrong orientation of the tip may lead to trauma of the urinary tract.

FEMALE

- With one hand spread the outer lips of the vulva to locate the urethral opening. With the other hand carefully introduce the catheter into the urethra until it reaches the bladder and the urine drains.

MALE / FEMALE

- Manual pressure on the lower abdomen towards the end of catheterisation will ensure that the bladder has been completely emptied.
- Immediately after the end of bladder emptying, withdraw the catheter.
- To stop the catheter from dripping, turn up the connector end of the catheter when removing.

Duration of use

For single use only.

Warning: Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination of the device may lead to infection. Catheter must be removed immediately after completion of urinary drainage.

Disposal

Discard the catheter according to the current local regulations.

Storage and handling conditions

Keep away from sunlight and heat. Keep dry. The product must be stored flat.

All Actreen® Hi-Lite Cath products are foldable, except CH18 and Actreen® Hi-Lite Cath 20cm.

Warning: Do not store the product (unit) in contact with sharp materials. Preferably use the discreet pouch for transport.

Warning: The expiry date only applies when Actreen® Hi-Lite Cath is correctly stored in an undamaged packaging.

Do not use after the expiry date.

Clinical benefits & Performances

Actreen® Hi-Lite Cath allows the drainage of urine from the bladder through the urethra.

Notice to the user

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

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DA B. Braun Medical A/S
FI B. Braun Medical Oy
NO B. Braun Medical A/S - Mathilde Henriksens vei 1
SE B. Braun Medical AB
B. Braun Avitum Italy S.p.A. | Via XXV Luglio 11 | 41037 Mirandola (MO) | ITALY