AT53102

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EN

Thank you for purchasing our product. Please read this user manual carefully before using the device. The medical device must not be used before reading and fully understanding these instructions. If you do not understand any warnings, cautions, or recommendations, contact a healthcare professional or the retailer to avoid product damage or personal injury. CAUTION: Check all parts for any damage caused during transport. If you notice any such damage, DO NOT USE the product. For more information, please contact the retailer.

INTENDED USER GROUP

This medical device is intended for persons suffering from diseases, dysfunctions, or injuries requiring treatment, rehabilitation, or compensation using this device (see the Intended Use section of this manual). The product may be purchased independently by the user or based on the recommendation of a doctor, therapist, or other specialist. In either case, available sizes, required functions, product variants, indications and contraindications, as well as the manufacturer's information, should be considered before use.

INTENDED USE

The forefoot offloading shoe is intended for temporary use to stabilize and protect the forefoot after surgical procedures.

Its construction transfers pressure to the heel, supporting the healing process of wounds, ulcers, and injuries in the front part of the foot.

The adjustable strap system allows the product to fit feet with dressings, bandages, or swelling.

INDICATIONS FOR USE

The orthosis is indicated for:

Fractures of the toes, Post-operative protection after surgeries in the forefoot area, Treatment of diabetic foot ulcers, particularly those located on the toes or forefoot, It can also be used for inflammatory conditions, injuries, or other pathologies that require pressure relief on the toes. Additionally, it supports the healing of wounds and skin lesions in the forefoot area.

CONTRAINDICATIONS

Physical or mental limitations (e.g., visual impairment) that prevent safe handling of the product.

Instability of the ankle joint requiring rigid stabilization of the entire foot and ankle.

Advanced foot deformities preventing proper fitting of the orthosis.

Inability to walk independently or with assistance while maintaining balance.

Severe ulcers or infections affecting the rear part of the foot or the heel.

LIST OF COMPONENTS: Forefoot offloading shoe, User manual

FEATURES

Material: 50% TPR, 30% polyester, 10% EVA, 10% nylon

Adjustable fit using two wide straps

Additional extension flaps to accommodate bandages, dressings, and swelling

Suitable for both left and right foot (neutral base)

Reduces painful pressure on the toes and forefoot

Transfers weight to the heel – ideal for toe fractures and diabetic foot ulcers

Lightweight, offloading sole construction

Rigid and durable rubber sole

Deep tread pattern for enhanced grip and stability on various surfaces

SIZE: XS: 35–37 S: 37–39 M: 40–42 L: 42–44 XL: 44–46

Note: If you are between sizes, choose the larger size if swelling or dressings are present.

INSTRUCTIONS FOR USE

Loosen the fastening straps. Carefully insert your foot into the orthosis, ensuring the heel is properly positioned. Adjust the straps for a snug but comfortable fit — the orthosis should be secure but not overly tight. Fasten the Velcro straps to ensure stable immobilization of the foot. If the top strap is too long, you may safely trim the excess.

WARNING – REPORTING OF A SERIOUS INCIDENT

In the event of a “serious incident” related to the use of this medical device that has led, could have led, or might lead to any of the following:

- a) death of the patient, user, or another person, or
- b) temporary or permanent deterioration of the health of the patient, user, or another person, or
- c) serious risk to public health,

the above serious incident must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides.

CAUTION: If you experience pain, allergic reactions, or any other unusual or unclear symptoms during the use of this medical device, contact a healthcare professional immediately. It is strictly forbidden to use the product in any manner other than its intended use.

CLEANING: Clean the surface of the shoe with a damp cloth and mild detergent (e.g., soap).

Allow to air dry completely after cleaning.

Do not machine wash or submerge the shoe in water. **STORAGE:** Store the product in a dry place away from direct sunlight.

DISPOSAL : After the product has reached the end of its service life, dispose of it as regular household waste. Always follow local waste separation regulations. Help protect the environment.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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