8.03.04.0122 - GS

Issue date: 2012-09-04

EN

Medical Sodium Hyaluronate Gel (for Bone Joint)

(Primary ingredients and performance)
SINGLOINT is a sterile gel solution containing disodium sodium hyaluronate or hyaluronic acid (hereinafter as "HA") in phosphate buffered saline, which is contained in a pre-filled syringe SINGLOINT for single intra-articular injection. The product is a unique type of single macromolecular hyaluronate, comprising repeating disaccharide units of glucosamine acid and N-acetyl glucosamine in chemical formula is (C₁₄H₂₁NO₁₁).

SINGLOINT is non-toxic animal, with the advantages of high viscosity, lubricity, physical absorbability and good biocompatibility. It is a medical product extracted and refined by high-tech biotechnology from animal-derived fermentation isolates. It is sterile, and free from pyrogenicity, allergy, genetic toxicity, or skin irritation. Due to its permeation protection of HA, the product has no risk of being a virus carrier and dissemination medium when containing animal tissues or organs. This product can be used as an analgesic agent to improve the joint function. When used in injections in joint cavity, due to its macromolecular characteristic and good viscoelasticity, it can effectively improve to produce high molecular chondrocyte extracellular matrix, which is beneficial for cartilage joint pain, increasing mobility, decreasing synovitis and delaying progression of disease.

(Indications)
SINGLOINT is indicated for the treatment of symptoms of osteoarthritis, joint pain. By replacing and supplementing the physiological function that is the cartilage joint, it reduces pain and improves joint function.

(Use)
A. Use: Attach needle to syringe.
1. Remove the top cap. Hold syringe and pull up cap of the syringe.
2. Insert needle. Hold the syringe body and firmly insert the hub of the needle/cannula a suitable penetration needed into the backside of the end of the syringe.
3. Tighten the needle. Tighten the needle by turning it firmly in a clockwise direction until it is seated in the proper position.

B. Injection instructions
1. As SINGLOINT is administered by intra-articular injection, the injection of SINGLOINT should be performed only by trained physician.
2. Aseptic injection procedure shall be followed.
3. After injection, the patient should be observed before performing the injection. The dosage regimen is to inject the product into the affected synovial joint cavity by one injection.
4. Discard the syringe and needle after the single use.

The dosage regimen is injecting into the affected synovial joint space once a week with a total of three injections.

(Adverse Reactions)
Mild transient local swelling, erythema and joint stiffness at the injection site are reported after intra-articular injections. These symptoms are generally mild and transient. They usually resolve within 72 hours without medical intervention. If symptoms are persistent, please consult with your treatment doctor.

(Precautions)
1. Strict aseptic technique must be followed for preventing patients from Cross infection.
2. The pre-filled syringe is intended for single use. The contents of the syringe must be used immediately after opening the packaging. Discard any rest of SINGLOINT. Secondary use is possible to local sterile infectious diseases.
3. Do not use if the product is in a damaged condition, or if any cracks or leakage appears on the pre-filled syringe.
4. Do not reuse SINGLOINT once used.
5. Do not inject SINGLOINT out of joint cavity or into the synovial tissue or capsule.
6. If joint infection present, infection should be eliminated before treating SINGLOINT.
7. Do not use SINGLOINT in patients with known hypersensitivity or allergy, especially to skin products, since sodium hyaluronate may precipitate when meet the ingredients.

(Contraindications)
1. SINGLOINT is contraindicated for patients with known hypersensitivity (allergy) to SINGLOINT or any other sodium hyaluronate.
2. SINGLOINT is contraindicated for patients with skin diseases or infection in the area of injection site.
3. SINGLOINT is contraindicated for patients with acute and chronic progressive arthritis, synovitis and joint tuberculosis.

(Dose)
The content for medical sodium hyaluronate is composed by three doses. Please use carefully. Do not use after the expiry date printed on the packaging. The expiry date refers to the product in its original packaging at the temperature of 2°C~20°C.
Keep out of reach of children.
The safety and effectiveness of SINGLOINT has been tested on pediatric patients, pregnant or lactating women.

(Storage)
Store between 2°C and 20°C, out of freeze.
Relative humidity: relative to no more than 80%.
Store in a clean, well-ventilated condition, avoid moisture contact.
Protect from light.

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