## Cutimed® Sorbact Pad®

(en)

Read the instructions for use before starting treatment with Cutimed® Sorbact Pad®.

### Device Description

Cutimed® Sorbact Pad® is a bacteria and fungi binding wound dressing, based on Sorbact® Technology. It consists of a Sorbact® DACC™-coated wound contact layer combined with an absorbent core. The product absorbs and retains exudate. Cutimed® Sorbact Pad® prevents and treats wound infections and facilitates the wound healing process. Sorbact® wound contact layer is shown to bind wound related pathogens up to seven days *in vitro*.

### Intended Purpose

Cutimed® Sorbact Pad® is intended for use in management of clean, contaminated, colonized or infected moderate to high exuding wounds, such as surgical wounds, traumatic wounds, pressure ulcers, diabetic foot ulcers and leg ulcers. Cutimed® Sorbact Pad® is intended to be used on superficial wounds.

### Intended User and Use Environment

Cutimed® Sorbact Pad® is intended to be used on children and adults. Cutimed® Sorbact Pad® is intended to be used by healthcare professionals or by lay persons under supervision of a healthcare professional. The dressing is intended to be used in healthcare facilities and home environment.

### Mode of Action

Cutimed® Sorbact Pad® binds common wound microorganisms, such as *Staphylococcus aureus* (including MRSA), *Streptococcus* spp., *Escherichia coli*, *Pseudomonas aeruginosa* and *Candida albicans*, as shown *in vitro*. There is no release of antimicrobial agents to the wound.

### Clinical Benefit

Use of Cutimed® Sorbact Pad® can reduce bioburden and prevent infection. The dressing helps reducing odor and pain, and can improve wound healing and decrease wound size.

### Contraindications

There are no known contraindications to the use of Cutimed® Sorbact Pad®. Warnings and precautions should be noted.

### Warnings and Precautions

Do not re-use. Cutimed® Sorbact Pad® is for single patient and single use only. Re-use may lead to cross contamination and increased risk of infection. Do not use if the peel pouch is already opened or damaged, since the sterility cannot then be guaranteed. Do not use on patients with a known hypersensitivity to the dressing (dialkylcarbamoyl chloride [DACC™], cellulose acetate, cellulose, cotton, polyester and viscose). Do not use in combination with fatty products, such as ointments, creams and solutions, as they may decrease the binding of microorganisms. Do not re-sterilize. If the treated condition deteriorates, fails to improve or a side effect is observed, consult an appropriate clinician.

### Instructions for Use

- Prepare the wound and surrounding skin according to local clinical practice.
- Select an appropriate dressing size for the wound. The wound pad should overlap the wound margins by at least 2 cm.
- Do not cut the dressing.
- Remove the dressing from the peel pouch using an aseptic technique.
- Apply the dressing. Ensure that the green contact layer comes into direct contact with the complete wound surface to allow microorganisms to bind to the dressing.
- Fixate with an elastic bandage or skin friendly tape.
- The dressing change frequency depends on exudate levels and the overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days.
- Should the dressing adhere to the wound, moisten the dressing to assist removal and to avoid disruption of the healing wound

Cutimed® Sorbact Pad® can be used in conjunction with compression therapy. Cutimed® Sorbact Pad® is MRI safe. Prior to radiation therapy, remove Cutimed® Sorbact Pad® if it is placed in the radiation field. A new dressing can be applied following treatment.

### Storage and Disposal

Cutimed® Sorbact Pad® shall be stored dry and kept away from sunlight. Disposal should be made according to local environmental procedures.

### Notice to User

Any serious incident that has occurred in relation to the device should be reported to ABIGO Medical AB (complaints@abigo.com) and the competent authority of your state.

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|                                                                                   |                                                                                    |                                       |                                                                                                                                                                                                                                                                                                                                                                                                                |
|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|---------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Artwork ID:<br>JBR19189.01                                                         | Review Loop:<br>01                    | <b>Printing Colours:</b><br> Cyan<br> Magenta<br> Yellow<br> Black |
|                                                                                   | Project Title:<br>C3132_GoPro Cutimed Sorbact Pad_IFU ROW                          | Date of Review Loop:<br>07.03.2024_BL |                                                                                                                                                                                                                                                                                                                                                                                                                |
|                                                                                   | Contact<br>ArtworkManagementMedical@essity.com                                     | Replaces Artwork ID:<br>JBR17462.02   |                                                                                                                                                                                                                                                                                                                                                                                                                |
|                                                                                   | Contact Address: BSN medical GmbH<br>Schützenstraße 1–3<br>22761 Hamburg   Germany | Technical Drawing for print ID<br>-   |                                                                                                                                                                                                                                                                                                                                                                                                                |
|                                                                                   | Languages:<br>en fr es it nl pt pt (br) el sl id ar cs sk ro                       | Note                                  | <b>Guiding Colours:</b><br> Diecut - DO NOT PRINT                                                                                                                                                                                                                                                                           |

## ESKO WebCenter Approval Information

**Project:** JBR19189.01  
**Project Number:** C3132 Cutimed Sorbact Pad  
**Project Name:** GoPro Cutimed Sorbact Pad  
**Reason for Change:** name change from CMED Sorbact dressing pad to CMED Sorbact Pad  
**Approval Email:** Micael Lagerqvist:micael.lagerqvist@essity.com, Inger Stenseke Larsson:inger.stenseke.larsson@essity.com, Pernilla Bottner:pernilla.bottner@essity.com  
**Approval Comment:** All below mentioned approvers attest to the accuracy and integrity of this document by signing the document with digital signature implemented in WebCenter.  
**Document:** JBR19189.01.pdf  
**Document description:**  
**Version:** 1  
**Version ID:** 00002\_0000035211  
**Downloaded:** Mar 19, 2024 at 7:42 AM GMT+01:00  
**By:** Heiko Hess (DEHHEEE)  
**Approval state:**

# Approved

| St<br>ag<br>e<br>Nr | Approver                                                                   | Role / Group                   | Date and<br>Time                         | Action                          | Comment |
|---------------------|----------------------------------------------------------------------------|--------------------------------|------------------------------------------|---------------------------------|---------|
| 1                   | Lagerqvist, Micael<br>(SEHMILG),<br>ESSITY - regular<br>user (Germany)     | PACKAGING<br>ROLE              | Mar 7, 2024 at<br>2:47 PM<br>GMT+01:00   | Approved<br>(Final<br>approval) |         |
| 1                   | Inger Stenseke<br>Larsson (SEHINST),<br>ESSITY - regular<br>user (Germany) | MARKETING<br>ROLE              | Mar 8, 2024 at<br>11:09 AM<br>GMT+01:00  | Approved<br>(Final<br>approval) |         |
| 1                   | Pernilla Bottner<br>(SEHPEBT),<br>ESSITY - regular<br>user (Germany)       | REGULATOR<br>Y AFFAIRS<br>ROLE | Mar 15, 2024<br>at 12:53 PM<br>GMT+01:00 | Approved<br>(Final<br>approval) |         |