Cutimed® Sorbion® Sorbact®

(EN)

Read the instructions for use before starting treatment with Cutimed® Sorbion® Sorbact®.

Device Description

Cutimed® Sorbion® Sorbact® is a bacteria and fungi binding wound dressing, based on Sorbact® Technology. It consists of a Sorbact® DACCTM-coated wound contact layer combined with a superabsorbent core and a water repellent backing, which prevents exudate strike through. The product absorbs and retains exudate, thereby reducing the risk of maceration. Cutimed® Sorbion® Sorbact® prevents and treats wound infections and facilitates the wound healing process. Sorbact® wound contact layer is shown to bind wound related pathogens up to seven days *in vitro*.

Intended Purpose

Cutimed® Sorbion® Sorbact® is intended for use in management of clean, contaminated, colonized or infected high to excessive exuding wounds, such as pressure ulcers, diabetic foot ulcers, leg ulcers, traumatic wounds and surgical wounds. Cutimed® Sorbion® Sorbact® is intended to be used on superficial wounds.

Intended User and Use Environment

Cutimed® Sorbion® Sorbact® is intended to be used on children and adults. Cutimed® Sorbion® Sorbact® is intended to be used by healthcare professionals or by lay persons under supervision of a healthcare professional. The dressing is intended to be used in healthcare facilities and home environment.

Mode of Action

Cutimed® Sorbion® Sorbact® binds common wound microorganisms, such as *Staphylococcus aureus* (including MRSA), *Streptococcus* spp., *Escherichia coli*, *Pseudomonas aeruginosa* and *Candida albicans*, as shown *in vitro*. There is no release of antimicrobial agents to the wound.

Clinical Benefit

Use of Cutimed® Sorbion® Sorbact® can reduce bioburden and prevent infection. The dressing helps reducing pain and can improve wound healing and decrease wound size. Cutimed® Sorbion® Sorbact® absorbs and retains high to excessive exudate levels.

Contraindications

There are no known contraindications to the use of Cutimed® Sorbion® Sorbact®. Warnings and precautions should be noted.

Warnings and Precautions

Do not re-use. Cutimed® Sorbion® Sorbact® is for single patient and single use only. Re-use may lead to cross contamination and increased risk of infection. Do not use if the peel pouch is already opened or damaged, since the sterility cannot then be guaranteed. Do not use on patients with a known hypersensitivity to the dressing (dialkylcarbamoyl chloride [DACCTM], cellulose acetate, copolyamide, polypropylene, cellulose, polyacrylate, synthetic rubber, polyurethane, acrylic print). Do not use in combination with fatty products, such as ointments, creams and solutions, as they may decrease the binding of microorganisms. Do not use on eyes, mucous membranes or in wound cavities. The dressing swells considerably after fluid absorption. Do not use on arterial bleeds or heavily bleeding wounds as the lost blood volume will be undetected. Do not re-sterilize. If the treated condition deteriorates, fails to improve or a side effect is observed, consult an appropriate clinician.

Instructions for Use

- Prepare the wound and surrounding skin according to local clinical practice.
- Select an appropriate dressing size for the wound. The wound pad should overlap the wound margins by at least 2cm. Do not cut the dressing.
- Remove the dressing from the peel pouch using an aseptic technique.
- Apply the dressing. Ensure that the green contact layer comes into direct contact with the complete wound surface to allow microorganisms to bind to the dressing.
- Fixate with an elastic bandage or skin friendly tape. Ensure the dressing is allowed to expand under the fixation.
- The dressing change frequency depends on exudate levels and the overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days.
- Should the dressing adhere to the wound, moisten the dressing to assist removal and to avoid disruption of the healing wound.

Cutimed® Sorbion® Sorbact® can be used in conjunction with compression therapy. Cutimed® Sorbion® Sorbact® is MRI safe. Prior to radiation therapy, remove Cutimed® Sorbion® Sorbact® if it is placed in the radiation field. A new dressing can be applied following treatment.

Storage and Disposal

Cutimed® Sorbion® Sorbact® shall be stored dry and kept away from sunlight. Store below 40°C/104°F. Disposal should be made according to local environmental procedures.

Notice to User

Any serious incident that has occurred in relation to the device should be reported to ABIGO Medical AB (complaints@abigo.com) and the competent authority of your state.

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		Artwork ID: JBR18331.01	Review Loop: 06	Printing Colours: Cyan Magenta Yellow Black
Project Title: C2628_Cutimed_Sorbion_Sorbact			Date of Review Loop: 30.01.2023 AF_KLND	
Contact ArtworkManagementMedical@essity.com			Replaces Artwork ID: JBR12204.04	
Contact Address: BSN medical GmbH Schützenstraße 1 – 3 22761 Hamburg Germany			Technical Drawing for print ID	
Languages: EN FR ES IT NL PT PT (Brazil) EL (GR) SL ID AR CZ SK			Note	

ESKO WebCenter Approval Information

Project: JBR18331.01
Project Name: Cutimed Sorbion Sorbact
Approval Email: Eva-Lie Essman:evalie.essman@essity.com, Jeanette Jochnick:jeanette.jochnick@essity.com, Håkan Jakobsson:hakan.jakobsson@essity.com, Ralf Timm:ralf.timm@essity.com, Birgit Struensee:birgit.struensee.ext@essity.com, Thu Nguyen:thu.nguyen@essity.com
Approval Comment: All below mentioned approvers attest to the accuracy and integrity of this document by signing the document with digital signature implemented in WebCenter.
Reason for Change: MDD to MDR
Project Number: C2628
Document: JBR18331.01.pdf
Document description:
Version: 2
Version ID: 00002_0000011887
Downloaded: Feb 2, 2023 at 1:11 PM GMT+01:00
By: Hilke Meyer (DEHHIME)
Approval state:

Approved

St ag e Nr	Approver	Role / Group	Date and Time	Action	Comment
1	Eva-Lie Essman (SEHEVES), ESSITY - regular user (Germany)	ABIGO REGULATORY AFFAIRS ROLE	Jan 30, 2023 at 3:38 PM GMT+01:00	Approved (Final approval)	
1	Jochnick, Jeanette (SEHJEJC), ESSITY - regular user (Germany)	ABIGO R&D ROLE	Jan 30, 2023 at 3:41 PM GMT+01:00	Approved (Final approval)	
1	Håkan Jakobsson (SEHHKJA), ESSITY - regular user (Germany)	ABIGO SUPPLY ROLE	Jan 30, 2023 at 3:55 PM GMT+01:00	Approved (Final approval)	
1	Ralf Timm (DEHRATI), ESSITY - regular user (Germany)	ABIGO PACKAGING ROLE	Jan 31, 2023 at 11:52 AM GMT+01:00	Approved (Final approval)	
1	Birgit Struensee (DEHBIST), ESSITY - regular user (Germany)	ABIGO MARKETING ROLE	Jan 31, 2023 at 12:41 PM GMT+01:00	Approved (Final approval)	
1	Thu Nguyen (SEHTHING), ESSITY - regular user (Germany)	ABIGO FINAL APPROVAL ROLE	Jan 31, 2023 at 1:32 PM GMT+01:00	Approved (Final approval)	