Cutimed® Sorbact Contact®

(en)

Read the instructions for use before starting treatment with Cutimed® Sorbact Contact®.

Device Description

Cutimed® Sorbact Contact® is a DACC™ -coated, non-absorbent bacteria and fungi binding wound contact layer, based on Sorbact® Technology, that allows passage of wound exudate into a secondary dressing. Cutimed® Sorbact Contact® prevents and treats wound infections and facilitates the wound healing process.

Intended Purpose

Cutimed® Sorbact Contact® is intended for use in management of clean, contaminated, colonized or infected exuding wounds, such as surgical wounds, traumatic wounds, pressure ulcers, diabetic foot ulcers and leg ulcers. Cutimed® Sorbact Contact® can be used on both superficial and deep wounds.

Intended User and Use Environment

Cutimed® Sorbact Contact® is intended to be used on children and adults. Cutimed® Sorbact Contact® is intended to be used by healthcare professionals, or by lay persons under supervision of a healthcare professional. The dressing is intended to be used in healthcare facilities and home environment.

Mode of Action

Cutimed® Sorbact Contact® binds common wound microorganisms, such as *Staphylococcus aureus* (including MRSA), *Streptococcus* spp., *Escherichia coli*, *Pseudomonas aeruginosa* and *Candida albicans*, as shown *in vitro*. There is no release of antimicrobial agents to the wound.

Clinical Benefit

Use of Cutimed® Sorbact Contact® can reduce bioburden and prevent infection. The dressing helps reducing odor and pain, and can improve wound healing and decrease wound size.

Contraindications

There are no known contraindications to the use of Cutimed® Sorbact Contact®. Warnings and precautions should be noted.

Warnings and Precautions

Do not re-use. Cutimed® Sorbact Contact® is for single patient and single use only. Re-use may lead to cross contamination and increased risk of infection. Do not use if the peel pouch is already opened or damaged, since the sterility cannot then be guaranteed. Do not use on patients with a known hypersensitivity to the dressing (dialkylcarbamoyl chloride [DACC™], cellulose acetate). Do not use in combination with fatty products, such as ointments, creams and solutions as they may decrease the binding of microorganisms. Do not wrap around extremities, as circulation may be compromised. Do not apply Cutimed® Sorbact Contact® in direct contact with bone and tendons, exposed blood vessels, organs or nerves. Do not re-sterilize. If the treated condition deteriorates, fails to improve or a side effect is observed, consult an appropriate clinician.

Instructions for Use

- Prepare the wound and surrounding skin according to local clinical practice.
- Select an appropriate dressing size for the wound. The dressing can be unfolded into suitable size and may overlap the wound margins if needed.
- Remove the dressing from the peel pouch using an aseptic technique.
- If the dressing is cut, use an aseptic technique. Do not tear the dressing. Discard any unused dressing.
- Apply the dressing. Ensure that the dressing comes into direct contact with the complete wound surface to allow microorganisms to bind to the dressing. In deep wounds, allow and fluffy up one or more dressings and fill the wound. Avoid dense packing. The dressing can also be used as a liner.
- Apply a secondary dressing appropriate for the exudate level such as a foam dressing, absorbent pad or gelling fiber dressing, and fixate.
- The dressing change frequency depends on exudate levels and the overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days.
- Should the dressing adhere to the wound, moisten the dressing to assist removal and to avoid disruption of the healing wound.

Cutimed® Sorbact Contact® can be used in conjunction with compression therapy. Cutimed® Sorbact Contact® can be used in conjunction with NPWT. Follow the Instructions for Use of the NPWT system. Cutimed® Sorbact Contact® is MRI safe. Prior to radiation therapy, remove Cutimed® Sorbact Contact® if placed in the radiation field. A new dressing can be applied following treatment.

Storage and Disposal

Cutimed® Sorbact Contact® shall be stored dry and kept away from sunlight.

Disposal should be made according to local environmental procedures.

Notice to User

Any serious incident that has occurred in relation to the device should be reported to ABIGO Medical AB (complaints@abigo.com) and the competent authority of your state.

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	Artwork ID: JBR19145.01	Review Loop: 03	Printing Colours:  Cyan  Magenta  Yellow  Black
Project Title: C3132_GoPro Cutimed Sorbact Contact_IFU ROW		Date of Review Loop: 07.02.2024_RB	
Contact ArtworkManagementMedical@essity.com		Replaces Artwork ID: JBR17385.02	
Contact Address: BSN medical GmbH Schützenstraße 1–3 22761 Hamburg Germany		Technical Drawing for print ID -	Guiding Colours:  Diecut - DO NOT PRINT
Languages: EN FR ES IT NL PT PT (Brazil) EL (GR) SL ID AR CS SK		Note	

ESKO WebCenter Approval Information

Project: JBR19145.01
Project Number: C3132 Cutimed Sorbact Contact
Project Name: GoPro Cutimed Sorbact Contact
Approval Email: Inger Stenseke Larsson:inger.stenseke.larsson@essity.com, Eva-Lie Essman:evalie.essman@essity.com, Micael Lagerqvist:micael.lagerqvist@essity.com
Reason for Change: name change from swab to Contact
Approval Comment: All below mentioned approvers attest to the accuracy and integrity of this document by signing the document with digital signature implemented in WebCenter.
Document: JBR19145.01.pdf
Document description:
Version: 3
Version ID: 00002_0000034474
Downloaded: Feb 16, 2024 at 10:56 AM GMT+01:00
By: Heiko Hess (DEHHEEE)
Approval state:

Approved

Stage Nr	Approver	Role / Group	Date and Time	Action	Comment
1	Inger Stenseke Larsson (SEHINST), ESSITY - regular user (Germany)	MARKETING ROLE	Feb 8, 2024 at 5:06 PM GMT+01:00	Approved (Final approval)	
1	Eva-Lie Essman (SEHEVES), ESSITY - regular user (Germany)	REGULATOR Y AFFAIRS ROLE	Feb 9, 2024 at 9:13 AM GMT+01:00	Approved (Final approval)	
1	Lagerqvist, Micael (SEHMILG), ESSITY - regular user (Germany)	PACKAGING ROLE	Feb 15, 2024 at 7:39 AM GMT+01:00	Approved (Final approval)	