PAVIS S.p.A.

Sede / Headquarter:

Via G. Rossini, 11/B - 21020 Buguggiate (VA) - Italy

Tel. (+39) 0332 45 54 63

Fax (+39) 0332 45 97 48

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 (UK)
1 - INDICATIONS (#131) • Distal ulnar or radial fractures, wrist sprains, flexor and extensor tendon injuries • Tendinitis, osteoarthritis, arthritis • Post-traumatic / post-surgical / post-plaster cast "joint rest" (#133) • Trauma/sprains of the wrist, injury of the flexor and extensor tendons • Tendinitis, osteoarthritis, arthritis, carpal tunnel • Post-traumatic / post-surgical / post-plaster cast "joint rest" (also at night)

2 - CONTRA-INDICATIONS AND SIDE-EFFECTS The chosen model must respond to the needs identified by the orthopaedic specialist/technician depending on the clinical diagnosis. Skin redness caused by rubbing may occur if it is not properly adjusted. As this is a NON INVASIVE item, consisting primarily of non-allergenic materials, no side effects of any type have ever been reported. Any severe incidents occurring in relation to the device must be reported to the manufacturer and the competent authority of the Member State where the user and/or patient resides.
--

3 - MAINTENANCE AND CARE • Before washing, remove the stays and close hook and loop parts     This device is non-biodegradable. Even though it is made of natural fabrics, after use, avoid release the product to the environment. AT THE END OF ITS USABLE LIFE HAS TO BE DISCHARGED ALONG WITH SOLID ORDINARY WASTE.

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